

Cofilin-1

Catalog No: #AP81964



Package Size: #AP81964-1 50ug #AP81964-2 100ug #AP81964-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Cofilin-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Cofilin, non-muscle isoform |
| Accession No. | Uniprot ID:Q5E9F7 |
| Uniprot | Q5E9F7 |
| GeneID | 534553; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CFL1

Note: This product is for in vitro research use only