

Collagen alpha-1(XV) chain

Catalog No: #AP81965



Package Size: #AP81965-1 50ug #AP81965-2 100ug #AP81965-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XV) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:O35206
Uniprot	O35206
GeneID	12819;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Col15a1

Note: This product is for in vitro research use only