

Metal transporter CNNM4

Catalog No: #AP81975



Package Size: #AP81975-1 50ug #AP81975-2 100ug #AP81975-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metal transporter CNNM4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ancient conserved domain-containing protein 4,mACDP4,Cyclin-M4,Acdp4,Kiaa1592
Accession No.	Uniprot ID:Q69ZF7
Uniprot	Q69ZF7
GeneID	94220;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cnnm4

Note: This product is for in vitro research use only