

Coronin-1A

Catalog No: #AP81977



Package Size: #AP81977-1 50ug #AP81977-2 100ug #AP81977-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Coronin-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Coronin-like protein A,Clipin-A,Coronin-like protein p57,Tryptophan aspartate-containing coat protein,TACO,Coro1
Accession No.	Uniprot ID:O89053
Uniprot	O89053
GeneID	12721;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Coro1a

Note: This product is for in vitro research use only