

Cannabinoid receptor 1

Catalog No: #AP81982



Package Size: #AP81982-1 50ug #AP81982-2 100ug #AP81982-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cannabinoid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CB-R,CANN6,CNR
Accession No.	Uniprot ID:P21554
Uniprot	P21554
GeneID	1268;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CNR1

Note: This product is for in vitro research use only