

Cofilin-1

Catalog No: #AP81986



Package Size: #AP81986-1 50ug #AP81986-2 100ug #AP81986-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cofilin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cofilin, non-muscle isoform
Accession No.	Uniprot ID:P18760
Uniprot	P18760
GeneID	12631;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cfl1

Note: This product is for in vitro research use only