

Cytochrome P450 1B1

Catalog No: #AP81988



Package Size: #AP81988-1 50ug #AP81988-2 100ug #AP81988-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 1B1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CYP1B1,Cytochrome P450CMEF,Cytochrome P450EF,Cyp1-b1
Accession No.	Uniprot ID:Q64429
Uniprot	Q64429
GeneID	13078;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp1b1

Note: This product is for in vitro research use only