

Collagen alpha-1(XVIII) chain

Catalog No: #AP81990



Package Size: #AP81990-1 50ug #AP81990-2 100ug #AP81990-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XVIII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:P39061
Uniprot	P39061
GeneID	12822;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Col18a1

Note: This product is for in vitro research use only