

Cytochrome c oxidase subunit 4 isoform 1, mitochondrial

Catalog No: #AP81992

Package Size: #AP81992-1 50ug #AP81992-2 100ug #AP81992-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cytochrome c oxidase polypeptide IV, Cytochrome c oxidase subunit IV isoform 1, COX IV-1, COX4
Accession No.	Uniprot ID:P00423
Uniprot	P00423
GeneID	281090;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: COX4I1

Note: This product is for in vitro research use only