

High affinity copper uptake protein 1

Catalog No: #AP81998



Package Size: #AP81998-1 50ug #AP81998-2 100ug #AP81998-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High affinity copper uptake protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Copper transporter 1,CTR1,Solute carrier family 31 member 1,CTR1
Accession No.	Uniprot ID:Q8WNR0
Uniprot	Q8WNR0
GeneID	397186;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC31A1

Note: This product is for in vitro research use only