

## Coronin-1A

Catalog No: #AP81999



Package Size: #AP81999-1 50ug #AP81999-2 100ug #AP81999-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Coronin-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Coronin-like protein A,Clipin-A,Coronin-like protein p57,Tryptophan aspartate-containing coat protein,TACO,CORO1
Accession No.	Uniprot ID:Q92176
Uniprot	Q92176
GeneID	282196;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CORO1A

Note: This product is for in vitro research use only