

C-X-C motif chemokine 14

Catalog No: #AP82005



Package Size: #AP82005-1 50ug #AP82005-2 100ug #AP82005-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chemokine BRAK,MIP-2G,Small-inducible cytokine B14,MIP2G,NJAC,SCYB14,PSEC0212,UNQ240/PRO273
Accession No.	Uniprot ID:O95715
Uniprot	O95715
GeneID	9547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CXCL14

Note: This product is for in vitro research use only