

C-X-C motif chemokine 16

Catalog No: #AP82006



Package Size: #AP82006-1 50ug #AP82006-2 100ug #AP82006-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-X-C motif chemokine 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Scavenger receptor for phosphatidylserine and oxidized low density lipoprotein, SR-PSOX, Small-inducible cytokine B16, Transmembrane chemokine CXCL16, SCYB16, SRPSOX, UNQ2759/PRO6714
Accession No.	Uniprot ID:Q9H2A7
Uniprot	Q9H2A7
GeneID	58191;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CXCL16

Note: This product is for in vitro research use only