

Cathelin-related antimicrobial peptide

Catalog No: #AP82007



Package Size: #AP82007-1 50ug #AP82007-2 100ug #AP82007-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cathelin-related antimicrobial peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cramp, Cathelin-like protein, CLP, Cnlp, Cramp
Accession No.	Uniprot ID:P51437
Uniprot	P51437
GeneID	12796;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Camp

Note: This product is for in vitro research use only