

## Gap junction alpha-1 protein

Catalog No: #AP82009



Package Size: #AP82009-1 50ug #AP82009-2 100ug #AP82009-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gap junction alpha-1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Connexin-43, Cx43, Gap junction 43 kDa heart protein, Cxn-43
Accession No.	Uniprot ID:P08050
Uniprot	P08050
GeneID	24392;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Gja1

Note: This product is for in vitro research use only