

Gap junction delta-2 protein

Catalog No: #AP82010

Package Size: #AP82010-1 50ug #AP82010-2 100ug #AP82010-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gap junction delta-2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Connexin-36,Cx36,Gap junction alpha-9 protein,Gja9
Accession No.	Uniprot ID:O54851
Uniprot	O54851
GeneID	14617;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gjd2

Note: This product is for in vitro research use only