

# Cystatin-C

Catalog No: #AP82012



Package Size: #AP82012-1 50ug #AP82012-2 100ug #AP82012-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cystatin-C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cystatin-3,Gamma-trace,Neuroendocrine basic polypeptide,Post-gamma-globulin
Accession No.	Uniprot ID:P01034
Uniprot	P01034
GeneID	1471;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CST3

Note: This product is for in vitro research use only