

## Cytochrome c oxidase subunit 5A, mitochondrial

Catalog No: #AP82015



Package Size: #AP82015-1 50ug #AP82015-2 100ug #AP82015-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytochrome c oxidase subunit 5A, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytochrome c oxidase polypeptide Va
Accession No.	Uniprot ID:P20674
Uniprot	P20674
GeneID	9377;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: COX5A

Note: This product is for in vitro research use only