

ATP-dependent Clp protease proteolytic subunit, mitochondrial

Catalog No: #AP82018

Orders: order@signalwayantibody.com

Package Size: #AP82018-1 50ug #AP82018-2 100ug #AP82018-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	ATP-dependent Clp protease proteolytic subunit, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Endopeptidase Clp
Accession No.	Uniprot ID:Q2KHU4
Uniprot	Q2KHU4
GeneID	535981;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CLPP

Note: This product is for in vitro research use only