

## VE-Cadherin Mouse Monoclonal Antibody(3G8)

Catalog No: #44264

Package Size: #44264-1 50ul #44264-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	VE-Cadherin Mouse Monoclonal Antibody(3G8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3G8
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous VE-Cadherin protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of VE-Cadherin
Target Name	VE-Cadherin
Other Names	Cadherin-5 (7B4 antigen) (Vascular endothelial cadherin) (VE-cadherin) (CD antigen CD144)
Accession No.	Swiss-Prot:P33151Gene ID:1003
Uniprot	P33151
GeneID	1003;
Calculated MW	80-115kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Rat Heart Tissue using VE-Cadherin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Heart Tissue using VE-Cadherin Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only