

## C-X-C motif chemokine 16

Catalog No: #AP82024



Package Size: #AP82024-1 50ug #AP82024-2 100ug #AP82024-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-X-C motif chemokine 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Scavenger receptor for phosphatidylserine and oxidized low density lipoprotein,SR-PSOX,Small-inducible cytokine B16,Transmembrane chemokine CXCL16,Srpsox
Accession No.	Uniprot ID:Q8BSU2
Uniprot	Q8BSU2
GeneID	66102;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cxcl16

Note: This product is for in vitro research use only