

Cystatin-B

Catalog No: #AP82026



Package Size: #AP82026-1 50ug #AP82026-2 100ug #AP82026-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystatin-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cystatin-beta,Liver thiol proteinase inhibitor,Stefin-B,Cst6
Accession No.	Uniprot ID:P01041
Uniprot	P01041
GeneID	25308;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cstb

Note: This product is for in vitro research use only