

## Beta-defensin 126

Catalog No: #AP82027



Package Size: #AP82027-1 50ug #AP82027-2 100ug #AP82027-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-defensin 126
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta-defensin 26,DEFB-26,Defensin, beta 126,Epididymal secretory protein 13.2,ESP13.2,HBD26,C20orf8,DEFB26
Accession No.	Uniprot ID:Q9BYW3
Uniprot	Q9BYW3
GeneID	81623;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:DEFB126

Note: This product is for in vitro research use only