

T-cell surface glycoprotein CD4

Catalog No: #AP82028



Package Size: #AP82028-1 50ug #AP82028-2 100ug #AP82028-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	T-cell differentiation antigen L3T4,T-cell surface antigen T4/Leu-3,CD4
Accession No.	Uniprot ID:P06332
Uniprot	P06332
GeneID	12504;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name: Cd4

Note: This product is for in vitro research use only