

CD63 antigen

Catalog No: #AP82029



Package Size: #AP82029-1 50ug #AP82029-2 100ug #AP82029-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD63 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD63
Accession No.	Uniprot ID:Q9XSK2
Uniprot	Q9XSK2
GeneID	404156;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD63

Note: This product is for in vitro research use only