

Deoxycytidine kinase

Catalog No: #AP82031



Package Size: #AP82031-1 50ug #AP82031-2 100ug #AP82031-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Deoxycytidine kinase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | dCK |
| Accession No. | Uniprot ID:P27707 |
| Uniprot | P27707 |
| GeneID | 1633; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DCK

Note: This product is for in vitro research use only