

Cytochrome P450 1B1

Catalog No: #AP82033



Package Size: #AP82033-1 50ug #AP82033-2 100ug #AP82033-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 1B1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CYP1B1
Accession No.	Uniprot ID:Q16678
Uniprot	Q16678
GeneID	1545;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CYP1B1

Note: This product is for in vitro research use only