

Cytochrome P450 1A2

Catalog No: #AP82038



Package Size: #AP82038-1 50ug #AP82038-2 100ug #AP82038-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 1A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CYP1A2,Cholesterol 25-hydroxylase,Cytochrome P-448,Cytochrome P-450d,Cytochrome P450-D,Cyp1a-2
Accession No.	Uniprot ID:P04799
Uniprot	P04799
GeneID	24297;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp1a2

Note: This product is for in vitro research use only