

## Tyrosine-protein kinase CSK

Catalog No: #AP82040



Package Size: #AP82040-1 50ug #AP82040-2 100ug #AP82040-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tyrosine-protein kinase CSK
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-Src kinase,Protein-tyrosine kinase MPK-2,p50CSK
Accession No.	Uniprot ID:P41241
Uniprot	P41241
GeneID	12988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Csk

Note: This product is for in vitro research use only