

CD276 antigen

Catalog No: #AP82041



Package Size: #AP82041-1 50ug #AP82041-2 100ug #AP82041-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD276 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	B7 homolog 3,B7-H3, Costimulatory molecule,CD276,B7h3
Accession No.	Uniprot ID:Q7TPB4
Uniprot	Q7TPB4
GeneID	315716;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd276

Note: This product is for in vitro research use only