

T-cell surface glycoprotein CD5

Catalog No: #AP82042



Package Size: #AP82042-1 50ug #AP82042-2 100ug #AP82042-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lymphocyte antigen 1, Ly-1, CD5, Ly-1
Accession No.	Uniprot ID:P13379
Uniprot	P13379
GeneID	12507;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd5

Note: This product is for in vitro research use only