

Complement decay-accelerating factor, GPI-anchored

Catalog No: #AP82046



Package Size: #AP82046-1 50ug #AP82046-2 100ug #AP82046-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement decay-accelerating factor, GPI-anchored
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DAF-GPI,CD55,Cd55a,Daf,Daf1
Accession No.	Uniprot ID:Q61475
Uniprot	Q61475
GeneID	13136;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd55

Note: This product is for in vitro research use only