

Collagen alpha-1(XV) chain

Catalog No: #AP82055



Package Size: #AP82055-1 50ug #AP82055-2 100ug #AP82055-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XV) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:P39059
Uniprot	P39059
GeneID	1306;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COL15A1

Note: This product is for in vitro research use only