

Disabled homolog 2

Catalog No: #AP82056

Package Size: #AP82056-1 50ug #AP82056-2 100ug #AP82056-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Disabled homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Adaptor molecule disabled-2,Differentially expressed in ovarian carcinoma 2,DOC-2,Mitogen-responsive phosphoprotein,Doc2
Accession No.	Uniprot ID:P98078
Uniprot	P98078
GeneID	13132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dab2

Note: This product is for in vitro research use only