

Cysteine-rich secretory protein LCCL domain-containing 2

Catalog No: #AP82057

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82057-1 50ug #AP82057-2 100ug #AP82057-3 1mg

Description

Product Name	Cysteine-rich secretory protein LCCL domain-containing 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cysteine-rich secretory protein 11,CRISP-11,LCCL domain-containing cysteine-rich secretory protein 2,CRISP11,LCRISP2,UNQ2914/PRO1156/PRO9783
Accession No.	Uniprot ID:Q9H0B8
Uniprot	Q9H0B8
GeneID	83716;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CRISPLD2

Note: This product is for in vitro research use only