

Cystatin-9

Catalog No: #AP82059



Package Size: #AP82059-1 50ug #AP82059-2 100ug #AP82059-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cystatin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Testatin,Cresp
Accession No.	Uniprot ID:Q9Z0H6
Uniprot	Q9Z0H6
GeneID	13013;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cst9

Note: This product is for in vitro research use only