

Cytochrome P450 1B1

Catalog No: #AP82060



Package Size: #AP82060-1 50ug #AP82060-2 100ug #AP82060-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Cytochrome P450 1B1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | CYP1B1,Cytochrome P450RAP |
| Accession No. | Uniprot ID:Q64678 |
| Uniprot | Q64678 |
| GeneID | 25426; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp1b1

Note: This product is for in vitro research use only