

## 60 kDa heat shock protein, mitochondrial

Catalog No: #AP82065



Package Size: #AP82065-1 50ug #AP82065-2 100ug #AP82065-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	60 kDa heat shock protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	60 kDa chaperonin,Chaperonin 60,CPN60,Heat shock protein 60,HSP-60,Mitochondrial matrix protein P1,HSP60
Accession No.	Uniprot ID:P31081
Uniprot	P31081
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HSPD1

Note: This product is for in vitro research use only