

Calcium and integrin-binding protein 1

Catalog No: #AP82066



Package Size: #AP82066-1 50ug #AP82066-2 100ug #AP82066-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcium and integrin-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CIB,Calmyrin,DNA-PKcs-interacting protein,Kinase-interacting protein,KIP,Cib,Kip,Prkdcip
Accession No.	Uniprot ID:Q9R010
Uniprot	Q9R010
GeneID	81823;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cib1

Note: This product is for in vitro research use only