

Cytochrome P450 2E1

Catalog No: #AP82067



Package Size: #AP82067-1 50ug #AP82067-2 100ug #AP82067-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2E1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	4-nitrophenol 2-hydroxylase,CYP11E1,Cytochrome P450-J,Cytochrome P450RLM6,Cyp2e,Cyp2e-1
Accession No.	Uniprot ID:P05182
Uniprot	P05182
GeneID	25086;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp2e1

Note: This product is for in vitro research use only