

Versican core protein

Catalog No: #AP82071



Package Size: #AP82071-1 50ug #AP82071-2 100ug #AP82071-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Versican core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Chondroitin sulfate proteoglycan core protein 2,Chondroitin sulfate proteoglycan 2,Glial hyaluronate-binding protein,GHAP,Large fibroblast proteoglycan,PG-M,CSPG2
Accession No.	Uniprot ID:P81282
Uniprot	P81282
GeneID	282662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VCAN

Note: This product is for in vitro research use only