

Pro-opiomelanocortin

Catalog No: #AP82074



Package Size: #AP82074-1 50ug #AP82074-2 100ug #AP82074-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pro-opiomelanocortin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	POMC,Corticotropin-lipotropin,Pomc1
Accession No.	Uniprot ID:P01193
Uniprot	P01193
GeneID	18976;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pomc

Note: This product is for in vitro research use only