

Cartilage oligomeric matrix protein

Catalog No: #AP82075



Package Size: #AP82075-1 50ug #AP82075-2 100ug #AP82075-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cartilage oligomeric matrix protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	COMP,Thrombospondin-5,TSP5
Accession No.	Uniprot ID:P49747
Uniprot	P49747
GeneID	1311;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COMP

Note: This product is for in vitro research use only