

10 kDa heat shock protein, mitochondrial

Catalog No: #AP82079



Package Size: #AP82079-1 50ug #AP82079-2 100ug #AP82079-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	10 kDa heat shock protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hsp10,10 kDa chaperonin,Chaperonin 10,CPN10
Accession No.	Uniprot ID:Q64433
Uniprot	Q64433
GeneID	15528;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspe1

Note: This product is for in vitro research use only