

Cytochrome b-245 heavy chain

Catalog No: #AP82081



Package Size: #AP82081-1 50ug #AP82081-2 100ug #AP82081-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome b-245 heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CGD91-phox,Cytochrome b(558) subunit beta,Cytochrome b558 subunit beta,Heme-binding membrane glycoprotein gp91phox,NADPH oxidase 2,Neutrophil cytochrome b 91 kDa polypeptide,Superoxide-generating NADPH oxidase heavy chain subunit,gp91-1,gp91-phox,p22 phagocyte B-cytochrome,NOX2
Accession No.	Uniprot ID:P04839
Uniprot	P04839
GeneID	1536;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CYBB

Note: This product is for in vitro research use only