

CREB3 regulatory factor

Catalog No: #AP82082



Package Size: #AP82082-1 50ug #AP82082-2 100ug #AP82082-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CREB3 regulatory factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Luman recruitment factor, LRF, C5orf41
Accession No.	Uniprot ID:Q8IUR6
Uniprot	Q8IUR6
GeneID	153222;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CREBRF

Note: This product is for in vitro research use only