

Complement C4-B

Catalog No: #AP82084



Package Size: #AP82084-1 50ug #AP82084-2 100ug #AP82084-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement C4-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Basic complement C4,C3 and PZP-like alpha-2-macroglobulin domain-containing protein 3,CO4,CPAMD3
Accession No.	Uniprot ID:P0C0L5
Uniprot	P0C0L5
GeneID	100293534;721;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C4B

Note: This product is for in vitro research use only