

D(1A) dopamine receptor

Catalog No: #AP82086



Package Size: #AP82086-1 50ug #AP82086-2 100ug #AP82086-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D(1A) dopamine receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Dopamine D1 receptor,D1AR
Accession No.	Uniprot ID:Q95136
Uniprot	Q95136
GeneID	281125;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DRD1

Note: This product is for in vitro research use only