

## Dual specificity protein phosphatase 1

Catalog No: #AP82087



Package Size: #AP82087-1 50ug #AP82087-2 100ug #AP82087-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dual specificity protein phosphatase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mitogen-activated protein kinase phosphatase 1,MAP kinase phosphatase 1,Protein-tyrosine phosphatase 3CH134,Protein-tyrosine phosphatase ERP,3ch134,Mkp1,Ptpn10,Ptpn16
Accession No.	Uniprot ID:P28563
Uniprot	P28563
GeneID	19252;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dusp1

Note: This product is for in vitro research use only