

Dual specificity protein phosphatase 9

Catalog No: #AP82088



Package Size: #AP82088-1 50ug #AP82088-2 100ug #AP82088-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity protein phosphatase 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mitogen-activated protein kinase phosphatase 4,MAP kinase phosphatase 4,MKP4
Accession No.	Uniprot ID:Q99956
Uniprot	Q99956
GeneID	1852;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DUSP9

Note: This product is for in vitro research use only